



Product Datasheet

CD5 Antibody (orb1826070)

Description	The CD5 antigen is a transmembrane glycoprotein expressed on the surface of practically all mature h...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant protein corresponding to the extracellular domain of human CD5 was used as the immunogen for the recombinant CD5 antibody.
Storage	Aliquot the recombinant CD5 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P06127
Hazard Information	This recombinant CD5 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



IHC staining of FFPE human tonsil tissue...