



Product Datasheet

MUC5AC Antibody (orb1826089)

Description	Mucins are a group of high molecular weight glycoproteins consisting of a mucin core protein and O-l...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human MUC5AC protein was used as the immunogen for the MUC5AC antibody.
Storage	Aliquot the MUC5 [°] C antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P98088
Hazard Information	This MUC5AC antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human stomach tissue...



SDS-PAGE analysis of purified, BSA-free ...