



Product Datasheet

Calcitonin Antibody / CALCA / CGRP (orb1826116)

Description	Calcitonin is a 32 amino acid polypeptide hormone that preserves skeletal integrity and reduces bloo...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of human Calcitonin protein (around aa 1-141) was used as the immunogen for the Calcitonin antibody.
Storage	Aliquot the Calcitonin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P01258
Hazard Information	This Calcitonin antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1:10-1:20 for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human thyroid tissu...



IHC staining of FFPE human medullary car...