



Product Datasheet

Tryptophan Hydroxylase Antibody / TPH1 (orb1826125)

Description	Phenylalanine hydroxylase (PAH), tyrosine hydroxylase (TH) and tryptophan hydroxylase (TPH) comprise...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human TPH1 protein was used as the immunogen for the Tryptophan Hydroxylase 1 antibody.
Storage	Aliquot the Tryptophan Hydroxylase 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2, kappa
Clonality	Monoclonal
Uniprot ID	P17752
Hazard Information	This Tryptophan Hydroxylase 1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human kidney cancer...



SDS-PAGE analysis of purified, BSA-free ...