



Product Datasheet

VMAT2 Antibody / SLC18A2 (orb1826198)


Description	Neurotransmission depends on the regulated exocytotic release of chemical transmitter molecules. Thi...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of human SLC18A2 protein was used as the immunogen for the SLC18A2 antibody.
Storage	Aliquot the SLC18A2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2, kappa
Clonality	Monoclonal
Uniprot ID	Q05940
Hazard Information	This SLC18A2 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human brain tissue ...



IHC staining of FFPE human prostate tiss...



SDS-PAGE analysis of purified, BSA-free ...