


## Product Datasheet

### Fas Antibody / CD95 (orb1826219)

<b>Description</b>	mAb rB-R18 specifically recognizes CD95, also known as Fas, a transmembrane glycoprotein with a MW o...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FACS, IF, IHC-Fr
<b>Immunogen</b>	Recombinant human CD95 protein was used as the immunogen for the CD95 antibody.
<b>Storage</b>	Aliquot the CD95 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	1 mg/ml in 1X PBS; BSA free, sodium azide free
<b>Isotype</b>	Mouse IgG2a, kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P25445</b>
<b>Hazard Information</b>	This CD95 antibody is available for research use only.
<b>Dilution Range</b>	Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Immunohistochemistry (Frozen): 1-2ug/ml for 30 minutes at RT
<b>Expiration Date</b>	12 months from date of receipt.



SDS-PAGE analysis of purified, BSA-free ...