


Product Datasheet

Fas Antibody / CD95 (orb1826221)

Description	mAb rB-R18 specifically recognizes CD95, also known as Fas, a transmembrane glycoprotein with a MW o...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IF, IHC-Fr
Immunogen	Recombinant human CD95 protein was used as the immunogen for the CD95 antibody.
Storage	Aliquot the CD95 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a, kappa
Clonality	Recombinant
Uniprot ID	P25445
Hazard Information	This CD95 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Immunohistochemistry (Frozen): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



SDS-PAGE analysis of purified, BSA-free ...