



Product Datasheet

Retinol Binding Protein 4 Antibody / RBP4 (orb1826233)

Description	Retinol (Vitamin A) is transported in the blood bound to its carrier protein, retinol-binding protei...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein (within amino acids 29-148) from the human protein was used as the immunogen for the RBP4 antibody.
Storage	Aliquot the RBP4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P02753
Hazard Information	This RBP4 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human liver tissue ...



SDS-PAGE analysis of purified, BSA-free ...