
Product Datasheet

Retinol Binding Protein 4 Antibody / RBP4 (orb1826236)

Description

Retinol (Vitamin A) is transported in the blood bound


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

IHC staining of FFPE human kidney tissue...

Immunogen

A recombinant partial protein (within amino acids 29-148) from the human protein was used as the immunogen for the Retinol Binding Protein 4 antibody.

Storage

Aliquot the Retinol Binding Protein 4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2c, kappa

Clonality

Monoclonal

Uniprot ID

P02753

Hazard Information


This Retinol Binding Protein 4 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.



SDS-PAGE analysis of purified, BSA-free ...



Analysis of a HuProt (TM) microarray con...