


Product Datasheet

Perforin Antibody / PRF1 (orb1826261)

Description	Perforin is a pore-forming protein that leads to osmotic lysis of the target cells and subsequently ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant human Perforin-1 protein fragment (within amino acids 413-552) was used as the immunogen for the Perforin-1 antibody.
Storage	Aliquot the Perforin-1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P14222
Hazard Information	This Perforin-1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human spleen tissue...