



Product Datasheet

GAD65 Antibody / GAD2 (orb1826307)

Description	This mAb recognizes a protein of 65kDa, which is identified as glutamic acid decarboxylase 2 (GDA2)....
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	A recombinant fragment of human GAD2/GAD65 protein (within amino acids 1-200) was used as the immunogen for the GAD65 antibody.
Storage	Aliquot the GAD65 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	Q05329
Hazard Information	This GAD65 antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating),Western blot: 1-2ug/ml,Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human thyroid tissu...



SDS-PAGE analysis of purified, BSA-free ...