



Product Datasheet

GATA3 Antibody (orb1826315)

Description	GATA-3 (GATA binding protein 3) is a member of the GATA family of transcription factors. This 50kDa ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human GATA-3 protein was used as the immunogen for the GATA-3 antibody.
Storage	Aliquot the GATA-3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P23771
Hazard Information	This GATA-3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate tiss...



IHC staining of FFPE human breast carino...