



Product Datasheet

GATA3 Antibody (orb1826320)


Description	GATA-3 (GATA binding protein 3) is a member of the GATA family of transcription factors. This 50kDa ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein (within amino acids 150-250) from the human protein was used as the immunogen for the GATA-3 antibody.
Storage	Aliquot the GATA-3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P23771
Hazard Information	This GATA-3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human bladder tissu...



IHC staining of FFPE human breast carino...



IHC staining of FFPE human lung carcinom...