



Product Datasheet

Neuropilin 1 Antibody / NRP1 (orb1826346)


Description	Neuropilin is a type I transmembrane receptor that has been implicated in aspects of axon growth and...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein (within amino acids 700-900) from the human protein was used as the immunogen for the Neuropilin-1 antibody.
Storage	Aliquot the Neuropilin-1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	O14786
Hazard Information	This Neuropilin-1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human kidney tissue...



SDS-PAGE analysis of purified, BSA-free ...



Analysis of a HuProt (TM) microarray con...