

## Product Datasheet

### Cynomolgus TIGIT Protein, hFc Tag (orb1826712)

<b>Description</b>	Recombinant Cynomolgus TIGIT protein with C-terminal human Fc tag
<b>Reactivity</b>	Monkey
<b>Conjugation</b>	Unconjugated
<b>Form/Appearance</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Tag</b>	C-Human Fc tag
<b>Note</b>	For research use only
<b>Purity</b>	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
<b>Source</b>	Mammalian
<b>MW</b>	The protein has a predicted molecular mass of 46.6 kDa after removal of the signal peptide. The apparent molecular mass of cTIGIT-hFc is approximately 35-55 kDa due to glycosylation.
<b>Uniprot ID</b>	<a href="#">A0A7N9I9V4</a>
<b>Expiration Date</b>	6 months from date of receipt.