
Product Datasheet

Human MUC16(13810-14451) Protein, His Tag (orb1826750)

Description

Recombinant human MUC16(13810-

10×His tag

Reactivity

Human

Conjugation

Unconjugated

Form/Appearance

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Tag

C-10×His tag

Note

For research use only

Purity

The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Source

Mammalian

MW

The protein has a predicted molecular mass of 73.2 kDa after removal of the signal peptide. The apparent molecular mass of MUC16(13810-14451)-His is approximately 100-130 kDa due to glycosylation.

Uniprot ID[Q8WXI7](#)**Expiration Date**

6 months from date of receipt.