

Product Datasheet

Human IGSF11 Protein, hFc Tag (orb1826833)

| | |
|------------------------|---|
| Description | Recombinant Human IGSF11 Protein with C-terminal human Fc tag |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Form/Appearance | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Storage | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Tag | C-Human Fc Tag |
| Note | For research use only |
| Purity | The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining. |
| Source | Mammalian |
| MW | The protein has a predicted molecular mass of 49.4 kDa after removal of the signal peptide. The apparent molecular mass of IGSF11-hFc is approximately 55-70 kDa due to glycosylation. |
| Uniprot ID | Q5DX21 |
| Expiration Date | 6 months from date of receipt. |