

## Product Datasheet

### FZD10 Antibody (orb1826903)

<b>Description</b>	Anti-FZD10(tabituximab biosimilar) mAb
<b>Species/Host</b>	Chimeric
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC
<b>Target</b>	FZD10
<b>Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Note</b>	For research use only
<b>Application notes</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Isotype</b>	IgG1(E239D, M241L)
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">Q9ULW2</a>
<b>Dilution Range</b>	ELISA 1:5000-10000; Flow Cyt 1:100
<b>Expiration Date</b>	12 months from date of receipt.