

Product Datasheet

Human Interleukin 1 β , IL-1 β ELISA Kit (orb1826974)

Description	Human Interleukin 1 β , IL-1 β ELISA Kit
Reactivity	Human
Range	2.344 pg/ml-150 pg/ml
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	P01584
Sensitivity	0.586 pg/ml
Expiration Date	Please enquire.