

Product Datasheet

Mouse Interleukin 17,IL-17 ELISA Kit (orb1826978)

Description	Mouse Interleukin 17,IL-17 ELISA Kit
Reactivity	Mouse
Range	10.938 pg/ml-700 pg/ml
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	Q62386
Sensitivity	2.735 pg/ml
Expiration Date	Please enquire.