
Product Datasheet

Vinculin VCL Antibody (Monoclonal, hVIN-1) (orb18279)

Description

Vinculin VCL Antibody (Monoclonal, hVIN-1)

Species/Host

Mouse

Reactivity

Gallus, Human, Monkey, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ICC, IF, IHC-Fr, WB

Immunogen

Human vinculin, purified from uterus.

Form/Appearance

Lyophilized

Concentration

Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.

Isotype

Mouse IgG1

Clonality


Monoclonal

MW


116615 MW

Uniprot ID
[P85972](#)
Expiration Date


12 months from date of receipt.



WB analysis of Vinculin using anti-Vinculin...



IF analysis of Vinculin using anti-Vinculin...



IHC analysis of Vinculin using anti-Vinculin...