



Product Datasheet

Human Gremlin / GREM1 Protein (orb18285)

Description	Human Gremlin, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Lys 25 - Asp 184 (A...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Gremlin
Preservatives	50 mM Tris, 150 mM NaCl, pH7.5 with 1 mM EDTA
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 44.2 kDa. The protein migrates as 48-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human Gremlin, Fc Tag (orb18285) is expressed from human 293 cells (HEK293). It contains AA Lys 25 - Asp 184 (Accession # NP_037504).
MW	44.2 kDa
Uniprot ID	O60565
NCBI	NP_037504
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Gremlin prote...



Immobilized Human BMP-2 Protein, premium...