

Product Datasheet

Dengue Virus NS1 antibody (orb183053)

| | |
|----------------------------|--|
| Description | Dengue Virus NS1 antibody |
| Species/Host | Rabbit |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, ICC, IF, IHC-Fr, IHC-P, WB |
| Immunogen | KLH conjugated synthetic peptide derived from human DENV-2 NS1 (801-900/3391 aa) |
| Target | DENV-2 NS1 |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 40 kDa |
| Uniprot ID | P30026 |
| Expiration Date | 12 months from date of receipt. |