


Product Datasheet

DAPK2 antibody (orb18331)


Description	Goat polyclonal antibody to DAPK2
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	DAPK2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: Immunohistochemistry: In paraffin embedded Human Breast shows staining of the luminal side of the cell membranes of epithelial cells in the ducts. Recommended concentration, 2-4ug/ml.WB: Approx 40-45kDa bands seen in mouse brain extracts. Recommended for use at 1-3 µg/ml.
Protein Sequence	KALHPRRRSSTS
Clonality	Polyclonal
MW	42.9; 25.0
Entrez	23604
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



orb18331 staining (2ug/ml) of mouse brai...



orb18331 (2.5ug/ml) staining of paraffin...



Immunohistochemical staining of Human Th...