



Product Datasheet

eIF2B epsilon (phospho-Ser539) antibody (orb183322)

Description	eIF2B epsilon (phospho-Ser539) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human EFHC1 around the phosphorylation site of Ser539 (PD(p-S)R)
Target	phospho-eIF2B epsilon (Ser539)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80 kDa
Uniprot ID	Q13144
Dilution Range	IHC-P=1:100-500, WB=1:500-2000, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ce...



Immunohistochemical staining of human li...