

Product Datasheet

IGF2BP2 siRNA (Mouse) (orb1834945)

Description	siRNA to inhibit IGF2BP2 expression using RNA interference
Reactivity	Mouse
Tested Applications	RNA-I
Form/Appearance	Lyophilized powder
Storage	Shipped at 4 °C. Store at -20 °C for one year.
Note	For research use only
Application notes	We recommends transfection with 10 nM - 100 nM siRNA 48 to 72 hours prior to cell lysis. Before resuspending, briefly centrifuge the tube to ensure the lyophilized siRNA is at the bottom of the tube. Resuspend the siRNA oligos to an appropriate concentration with DEPC water. For example, resuspend one tube of 5 nmol siRNA oligo in 250 µl of DEPC water to get a final concentration of 20 µM.
Purity	> 97%
Source	Synthetic
Uniprot ID	Q5SF07
Entrez	319765
Expiration Date	12 months from date of receipt.