

Product Datasheet

TNFAIP8L3 siRNA (Mouse) (orb1835529)

Description	siRNA to inhibit TNFAIP8L3 expression using RNA interference
Reactivity	Mouse
Tested Applications	RNA-I
Form/Appearance	Lyophilized powder
Storage	Shipped at 4 °C. Store at -20 °C for one year.
Note	For research use only
Application notes	We recommends transfection with 10 nM - 100 nM siRNA 48 to 72 hours prior to cell lysis. Before resuspending, briefly centrifuge the tube to ensure the lyophilized siRNA is at the bottom of the tube. Resuspend the siRNA oligos to an appropriate concentration with DEPC water. For example, resuspend one tube of 5 nmol siRNA oligo in 250 µl of DEPC water to get a final concentration of 20 µM.
Purity	> 97%
Source	Synthetic
Uniprot ID	Q3TBL6
Entrez	244882
Expiration Date	12 months from date of receipt.