


Product Datasheet


FCRL4 antibody (orb183672)

Description	FCRL4 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic peptide (please contact us for details)
Target	FCRL4
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	57 kDa
Uniprot ID	P29351 , P41499
NCBI	031282 , 112572
Entrez	83417
Dilution Range	IHC-P: 1: 100-600, WB: 1:100-800
Expiration Date	12 months from date of receipt.




IHC-P image of pig stomach tissue using

...



Immunohistochemical staining of pig smal...



IHC-P image of mouse liver tissue using ...