



Product Datasheet

CNOT2 antibody (orb1837)

| | |
|----------------------------|---|
| Description | CNOT2 antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P, WB |
| Immunogen | KLH conjugated synthetic peptide derived from human CNOT2 (441-540/540 aa) |
| Target | CNOT2 |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 59 kDa |
| Uniprot ID | Q9NZN8 |
| Dilution Range | WB=1:500-2000, IHC-P=1:100-500 |
| Expiration Date | 12 months from date of receipt. |



IHC-P analysis of mouse bladder using CN...



Western blot: Lane 1:U87MG (Human) Cell;L...