



Product Datasheet

FOXP3 antibody (orb18370)


Description	Goat polyclonal antibody to FOXP3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC, WB
Target	FOXP3 / SCURFIN
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.IHC: In paraffin embedded Human Tonsil, stains nuclei of some cells in the interfollicular area. Recommended concentration, 0.3-1 ug/ml.
Protein Sequence	SQRPSRCSNPTPGP
Clonality	Polyclonal
MW	47.2; 43.4
Entrez	50943
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



orb18370 (1µg/ml) staining of paraffin ...



orb18370 staining of CD25-sorted (Treg) ...



Western blot analysis of Human Muscle(la...