



Product Datasheet

RAB11A antibody (orb18403)

Description	Goat polyclonal antibody to RAB11A
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	RAB11A / YL8
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 26kDa band observed in Mouse Heart lysates (calculated MW of 24.4kDa according to NP_004654.1). In transfected HEK293 transiently expressing RAB11A a band of approx. 28kDa is observed. This band is not observed in the non-transfected HEK293. Recommended for use at 1-3 µg/ml.
Protein Sequence	TENKPKVQCCQNI
Clonality	Polyclonal
MW	24.4
Entrez	8766
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb18403 staining (2ug/ml) of mouse hear...



HEK293 overexpressing RAB11A(RC200352) a...