

Product Datasheet

Histone H3 (phospho-Thr45) antibody (orb184169)

Description	Histone H3 (phospho-Thr45) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Histone H3 around the phosphorylation site of Thr45 (41-136/136 aa)
Target	phospho-Histone H3 (Thr45)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17 kDa
Uniprot ID	Q16695
Expiration Date	12 months from date of receipt.