

Product Datasheet

Histone H4 (phospho-Ser1) antibody (orb184172)

Description	Histone H4 (phospho-Ser1) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Histone H4 around the phosphorylation site of Ser1 ((p-S)GR)
Target	phospho-Histone H4 (Ser1)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11 kDa
Uniprot ID	P62805
Expiration Date	12 months from date of receipt.