



Product Datasheet

His Tag antibody (orb184194)

Description	His Tag antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated HHHHHH peptide
Target	His Tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:2000-5000
Expiration Date	12 months from date of receipt.



Western blot: His-tagged FMDV VP1 protein...



Western blot: Lane 1: His-Tagged Fusion Pr...