
Product Datasheet

IL13 receptor alpha 1 (phospho-Tyr405) antibody (orb184426)

Description

IL13 receptor alpha 1 (phospho-

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested Applications

WB

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human IL13 receptor alpha 1 around the phosphorylation site of Tyr405 (DI(p-Y)EK)

Target

phospho-IL13RA1 (Tyr405)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

42 kDa

Uniprot ID[P78552](#)**Dilution Range**

WB=1:500-2000

Expiration Date

12 months from date of receipt.