

Product Datasheet

KAT13C (phospho-SER736) antibody (orb184812)

Description	KAT13C (phospho-SER736) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human KAT13C around the phosphorylation site of Ser736 (PV(p-S)PK)
Target	phospho-KAT13C (SER736)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	159 kDa
Uniprot ID	Q15596
Expiration Date	12 months from date of receipt.