

Product Datasheet

UBE2T antibody (orb18482)

Description	Goat polyclonal antibody to UBE2T
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	HSPC150 / UBE2T
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx. 22kDa band observed in lysates of cell lines K562, HEK293 and A549 (calculated MW of 22.5kDa according to NP_054895.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	QLVGIEKKFHPDV
Clonality	Polyclonal
MW	22.5
Entrez	29089
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb18482
(0.1µg/ml)
staining of
K562 lys...