

Product Datasheet

MNK1 (phospho-Thr385) antibody (orb185175)

Description	MNK1 (phospho-Thr385) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MNK1 around the phosphorylation site of Thr385 (LP(p-T)P)
Target	phospho-MNK1 (Thr385)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51 kDa
Uniprot ID	Q9BUB5
Expiration Date	12 months from date of receipt.