

Product Datasheet

TCF19 antibody (orb18576)

Description	Goat polyclonal antibody to TCF19
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TCF19 / SC1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 35 kDa band observed in 293 lysate (predicted MW of 41kDa according to NP_009040). Recommended for use at 1-2 µg/ml.
Protein Sequence	RSTAKAPSDTPAHE
Clonality	Polyclonal
MW	38.3
Entrez	6941
Dilution Range	ELISA: 1:32000, WB: 1-2 µg/ml
Expiration Date	12 months from date of receipt.



orb18576 staining (1.5ug/ml) of 293 lysa...