

Product Datasheet

COP1 (phospho-Ser387) antibody (FITC) (orb186976)

Description	Rabbit polyclonal antibody against phospho-COP1 (Ser387) conjugated to FITC
Species/Host	Rabbit
Reactivity	Rat
Conjugation	FITC
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human COP1 around the phosphorylation site of Ser387 TA(p-S)QL
Target	COP1
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80kDa
Uniprot ID	Q8NHY2
Dilution Range	ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.