


Product Datasheet

ZMYND11 antibody (orb18744)

Description	Goat polyclonal antibody to ZMYND11
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	BS69 / ZMYND11
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Preliminary experiments gave no signal but low background in human heart and brain lysates at up to 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	SRVHGMHPKETT
Clonality	Polyclonal
MW	71.0; 66.5; 64.4; 60.9; 64.4; 60.1; 66.6
Entrez	10771
Dilution Range	ELISA: 1:32000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...