



Product Datasheet

ACTL6A / ACTL6B antibody (orb18766)


Description	Goat polyclonal antibody to ACTL6A / ACTL6B
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	BAF53A and BAF53B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000. WB: Approx 48kDa band observed in Hela lysates (predicted MW of 49kDa according to NP_057272, NP_004292). Recommended for use at 1-3 µg/ml.
Protein Sequence	YEEGGKQCVERKCP
Clonality	Polyclonal
MW	47.5; 43.2; 46.9
Entrez	51412 , 86
Dilution Range	ELISA: 1:16000, WB: 0.5-2µg/ml, IF/ICC: 10µg/ml, FACS: 10ug/ml
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of paraformalde...



Immunofluorescence analysis of paraforma...



Immunofluorescence analysis of paraforma...