


## Product Datasheet

# Goat Anti-Rabbit IgG (H&L)-HRP Polymer (Ready to use) (orb1878755)

<b>Description</b>	Ready to use Goat Polyclonal Secondary Antibody to Rabbit IgG (H&L) HRP polymer labeled
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, IH, WB
<b>Immunogen</b>	Rabbit IgG
<b>Form/Appearance</b>	Liquid in PBS, pH 7.2, containing 1% rAlbumin, 0.05% Proclin
<b>Storage</b>	Shipped and store at 4°C for one year.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Goat
<b>Dilution Range</b>	Ready to use
<b>Expiration Date</b>	12 months from date of receipt.



Immunohistochemical analysis staining in...