


Product Datasheet

Goat Anti-Mouse IgG (H&L)-HRP Polymer (Ready to use) (orb1878756)

Description	Ready to use Goat Polyclonal Secondary Antibody to Mouse IgG (H&L) HRP polymer labeled
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA, IH, WB
Immunogen	Mouse IgG
Form/Appearance	Liquid in PBS, pH 7.2, containing 1% rAlbumin, 0.05% Proclin
Storage	Shipped and store at 4°C for one year.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	Ready to use
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis staining in...