


Product Datasheet

Streptavidin-FITC labeled Antibody (orb1878761)

Description	Streptavidin to FITC labeled
Conjugation	FITC
Tested Applications	ELISA, FC, IF, IH
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.02% Sodium Azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Source	Native
Dilution Range	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)
Expiration Date	12 months from date of receipt.



Line colors represent the approximate vi...