

Product Datasheet

Rabbit Anti-Human IgG F(ab')₂ - APC (orb1878770)

Description	Rabbit Polyclonal Secondary Antibody to Human IgG F(ab') ₂ APC labeled
Reactivity	Human
Conjugation	APC
Tested Applications	ELISA, FC, IF, IH
Immunogen	Human IgG F(ab') ₂
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Dilution Range	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)
Expiration Date	12 months from date of receipt.