

## Product Datasheet

### Goat Anti-Rabbit IgG (H&L) - PE (orb1878853)

|                            |                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Goat Polyclonal Secondary Antibody to Rabbit IgG (H&L) PE labeled                                                    |
| <b>Reactivity</b>          | Rabbit                                                                                                               |
| <b>Conjugation</b>         | PE                                                                                                                   |
| <b>Tested Applications</b> | ELISA, FC, IF, IH                                                                                                    |
| <b>Immunogen</b>           | Rabbit IgG                                                                                                           |
| <b>Form/Appearance</b>     | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                     |
| <b>Note</b>                | For research use only                                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                                                           |
| <b>Source</b>              | Goat                                                                                                                 |
| <b>Dilution Range</b>      | IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                      |