

Product Datasheet

Goat Anti-Mouse IgM mu chain - PE (orb1878855)

Description	Goat Polyclonal Secondary Antibody to Mouse IgM mu chain PE labeled
Reactivity	Mouse
Conjugation	PE
Tested Applications	ELISA, FC, IF, IH
Immunogen	Mouse IgM mu chain
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)
Expiration Date	12 months from date of receipt.