

## Product Datasheet

### Goat Anti-Rat IgG (H&L) - Cy5.5 (orb1878865)

<b>Description</b>	Goat Polyclonal Secondary Antibody to Rat IgG (H&L) Cy5.5 labeled
<b>Reactivity</b>	Rat
<b>Conjugation</b>	Cy5.5
<b>Tested Applications</b>	ELISA, FC, IF, IH
<b>Immunogen</b>	Rat IgG
<b>Form/Appearance</b>	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Goat
<b>Dilution Range</b>	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)
<b>Expiration Date</b>	12 months from date of receipt.