

Product Datasheet

Goat Anti-Rat IgG (H&L) - Cy5 (orb1878866)

Description	Goat Polyclonal Secondary Antibody to Rat IgG (H&L) Cy5 labeled
Reactivity	Rat
Conjugation	Cy5
Tested Applications	ELISA, FC, IF, IH
Immunogen	Rat IgG
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)
Expiration Date	12 months from date of receipt.